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Caries Detector

As dentists, you can sometimes encounter difficulties in identifying caries. It is therefore reassuring to have access to a clinical aid for effective removal of carious dentine, such as Caries Detector.

Caries Detector is the excellent aid in conservative dentistry. It aids in the excavation of the outer decalcified and infected carious dentine layer, permitting an optimal caries removal.

It also helps to minimize the removal of remineralizable, healthy dentine, protecting the vitality of the pulp and as a result to conserve the maximal amount of healthy tissue.

Carious dentin consists of two layers. The outer layer is decalcified and infected with bacteria. This layer cannot be remineralized and must be removed. The inner layer, however, is only partially decalcified and is not infected. This layer is still remineralizable and should be preserved as much as possible.

Composition: Glycol- Dye

Intended use:

- ▶ Detecting carious dentine for cavity preparation

Advantages:

- ▶ optimal removal of the infected carious dentine;
- ▶ optimal conservation of healthy, remineralizable dentine, thus minimizing accidental pulpal injury

Composition:

- ▶ Isolate carious tooth with rubber dam. After initial excavation of carious enamel and dentine over infected site, wash and dry.
- ▶ Place one drop of CARIES DETECTOR on a sponge pledget and apply to the cavity.
- ▶ Leave for 10 seconds, then rinse with water. Remove the exposed outer carious dentin, which will be stained red.
- ▶ Repeat detection to ensure removal of the infected carious

Precautions:

- ▶ Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- ▶ If the product comes in contact with the human body, take the following actions:
 - ▶ If the product gets in the eye: Immediately wash the eye with copious amounts of water and consult a Doctor.
 - ▶ If the product comes in contact with the skin: Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
 - ▶ If the product comes in contact with the mucous membrane: Immediately wipe it off with cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
- ▶ Use caution when using the product to prevent swallowing.
- ▶ Do not use the same disposable brush for different patients to prevent cross contamination. Discard the brush after use.
- ▶ Avoid direct contact with the skin and / or soft tissue to prevent hypersensitivity.
- ▶ Wear gloves or take appropriate precautions when using the product.
- ▶ Do not apply the product directly from the bottle to carious cavity.

Note:

- ▶ Recommended method for removal of carious dentin: Only the non-sensitive stainable tissue should be removed, carefully avoiding cutting the sensitive unstainable tissue.
- ▶ A low-speed rotary instrument should be used for carious dentin removal since the friction heat generated by a high-speed tool may conduct pain sensation even before reaching the sensitive unstained tissue.

Cavity preparation procedure using the product:

- ▶ The product can be used in two ways:
 - ▶ The product is applied to the cavity area before opening the cavity. This procedure is suggested for a shallow open cavity because earlier disclosure of undermined free enamel may facilitate opening the cavity.
 - ▶ The product is used to disclose the residual outer carious dentin after opening a cavity and removing heavily discolored and highly softened dentin. This procedure is suggested for a deep cavity with a large amount of carious dentin for expedience. Careless gross reduction before application of the product, however, may cause pain from unintentional over-reduction.


Handling and manipulation precautions:

- ▶ The use of this product is restricted to a licensed dental professional.
- ▶ Provide moisture and contamination control with the use of a rubber dam.
- ▶ Do not mix the product with any other dental materials.
- ▶ Do not use the product for any purposes other than those specified in INDICATIONS.
- ▶ The product must be used by the expiration date indicated on the package.
- ▶ The product must be stored at 5-25°C when not in use.
- ▶ Keep away from extreme heat or direct sunlight

Contraindication:

- ▶ This product has no contraindications. Do not use only for people who are allergic to the raw materials of this product.

Reg. NO. at Medical equipment department C01855

 Consult instructions for use	 Indicates the batch code	 Catalogue number	 Indicates the date after which the medical device must not be used	 Indicates a medical device that needs to be protected from light sources	 Indicates a medical device that needs to be protected from moisture.
 Do not use if Package is Damaged	 Indicates the medical device (Name & Address)	 Non-sterile Indicates a non-sterile device (Name & Address)	 Indicates the device intended for use by a professional person only	 Indicates that the medical device must be protected from the use of steam	 5°C 25°C Indicates the Temperature Limit

ORVABON
MORVABON Co.
Produced by Morvabon Co.
www.morvabon.com
info@morvabon.com
Made in IRAN

Caries Detector

- ISO13485 : 2016 , US cert, IAS
- Iso 9001: 2015, Tuv certification
- Iso 14001: 2015, Tuv certification
- Iso 18001: 2007, Tuv certification
- Iso 10002: 2014, Tuv certification
- Iso 10004: 2012, Tuv certification

**MEDICAL
EQUIPMENT
DEPARTMENT**
Reg. NO. at Medical
equipment department
C01855

manufacturer
Pardis Morvabon Pharmed Co.:
Head Office: Unit 154, 7nd floor, 7 block,
Bojzaz Bd, Azadi Street, Tehran, Province, Iran
Tel: 0098 21 6658 1438
Fax: 0098 21 6690 8819
Factory: Unit 135, 136, 139, 2nd floor, Seraj Bd.,
20thkm of Damavand Rd., Tehran, Iran
Tel: 0098 21 7625 1702
Fax: 0098 21 7648 1457