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Liquid for root canals rinsing

Composition:

chlorhexidine Gluconate 2%, water Ph. Eur, Stabilizer

Intended use:

- ▶ Preparation of root canals
- ▶ Liquid is used for rinsing of root canals in order to remove debris
- ▶ In case of infection caused by Enterococcus faecalis or Candida albicans
- ▶ In case of non-healing inflammation periapical of tissues

Advantages:

- ▶ Longest and most effective antibacterial activity
- ▶ Elimination of organic compound and bacteria into the root canal
- ▶ Increase bond strength with dentin bonding agents
- ▶ More effective on E. faecalis than the Hypochlorite sodium %5.25 at both time periods
- ▶ Antimicrobial effect on staphylococcus aureus and candida albicans in 15 seconds
- ▶ The positive-charged molecules of CHX can adsorb onto dentine and prevent microbial colonization on the dentine surface for some time beyond the actual medication period
- ▶ Remains antimicrobial effect for least 72 hours and 2 weeks after it is placed in the canal

Precautions:

- ▶ Carefully read and understand instructions before using.
- ▶ It's necessary to use rubber-dam while working with Morva sept.
- ▶ Use eye protection, mask and gloves when using Morva sept
- ▶ Avoid contact with oral soft tissue, eyes or skin
- ▶ If accidental contact with soft tissues, eyes or skin occurs,
- ▶ immediately flush with plenty of water
- ▶ In case of consumption drink copious amount of water.
- ▶ Do not induce vomiting. Seek medical advice
- ▶ Do not inject this product

- ▶ For dental use only.
- ▶ At the end of the work, it is not necessary to wash with water or serum, and just remove the extra amounts from the patient's mouth.
- ▶ Do not mix directly with a NaOCl solution in the canal as a harmful brown precipitate will form.

warning:

- ▶ Do not use after expiration date
- ▶ Keep out of reach of children
- ▶ Protect from direct sunlight
- ▶ Store in the original container at room temperature


Direction for use:

- ▶ MORVA-Sept %2 solution is recommended for endodontic disinfection as final rinse and it should be used after smear layer removal by EDTA solution. This allows the CHX to come in contact with the bacteria on the walls of the canal and in the dentine tubules.
- ▶ Insert the liquid into the root canal using a syringe with a side vent needle during mechanical preparation. Rinse the canals copiously, approximately 10 ml for one canal. While rinsing root canals, it is necessary to rinse with physiological saline between rinsing sodium hypochlorite and chlorhexidine Gluconate. It prevents from precipitating of a tawny deposit.

Contraindication:

- ▶ This product is not contraindicated for children, pregnant women, the elderly and other specific diseases and can be used for all people who are not allergic to the active substance

Available packaging: Bottle of 250 ml and 1 lit

MORVA-Sept %2 Solution

Chlorhexidine Gluconate %2

- ISO13485 : 2016 , US cert, IAS
- Iso 9001: 2015, Tuv certification
- Iso 14001: 2015, Tuv certification
- Iso 18001: 2007, Tuv certification
- Iso 10002: 2014, Tuv certification
- Iso 10004: 2012, Tuv certification

Production
License No.
77340703
IMD: 70803369

manufacturer
Pardis Morvabon Pharmed Co.:
Head Office: Unit 154, 7th floor, 7 block,
Borjshaz Bd, Azadi Street, Tehran, Province, Iran
Tel: 0098 21 6658 1438
Fax: 0098 21 6690 8819
Factory: Unit 135, 136, 139, 2nd floor, Seraj Bd.,
20thkm of Damavand Rd.,, Tehran, Iran
Tel: 0098 21 7625 1702
Fax: 0098 21 7648 1457